


Rocket® Haemorrhoid Injection Sets




Rocket® Haemorrhoid Injection Sets provide a safe, reliable method of injection of oily Phenol BP or other sclerosing agents as an outpatient treatment for haemorrhoids.



The design features:

- **COMFORTABLE** three finger ring grip for ease of injection of viscid Phenol.
- **ANGLED or STRAIGHT** luer lock needle for ease security and improved ease of access
- **WIDE BORE** needle with 18G distal tip to provide maximum control during injection.
- **SUPPLIED COMPLETE** with filling cannula for ease of use in the busy clinic setting.
- **RELIABLE**, practical and competitively priced to ensure best value for your department.

Ordering Information		
Rocket® Haemorrhoid Injection Set, Supplied with 10ml, 3 ring syringe, haemorrhoid injection needle, filling cannula. Supplied sterile, for single use. 25 sets per carton	Straight, 10cm, 4"	856210
	Straight, 20cm 8"	859503
	Angled, 10cm, 4"	856211



North Island: 0800 658 814 South Island: 0800 804 546
 customerservices@uslmedical.co.nz