

# Rocket® HSG & Sonography Procedure Packs



Rocket HSG & Sonography Procedure Packs provide all the necessary single use items for cost effective uterine canalisation for hysterosalpingography and hysterosalpingo-contrast sonography.

## CONTENTS:

- Rocket® HSG Catheter
- Vaginal speculum with screw retention
- Sponge stick for prep solution.
- Sponge holder with 5 x Green, X-ray detectable swabs
- Tenaculum
- Cervical dilator
- Water soluble lubricant
- 4 x 30cm swab sticks
- Filling kwill
- 20ml & 3ml syringes
- 30cm contrast extension tube
- Sanitary towel

## Rocket® HSG Catheters have these well-established and trusted features:

- **LATEX FREE** balloon minimises risk of patient tissue reactions.
- **SMOOTH TIP & LOW PROFILE BALLOON** ensures low trauma insertion and removal – ideal for the nulliparous cervix.
- **REMOVABLE STIFFENING SHEATH** to aid insertion in the heavily retroverted or antverted cervix.
- **CONTRAST PATH TAP** ensures retention of contrast when further volumes of media are required.
- **INTEGRAL WIRE STYLET** for improved malleability without using tiresome removable stylets (ST versions)
  - ◇ Improves ease of insertion
  - ◇ Increased value for money
- **COMPATIBLE** with all aqueous based contrast media.
- **COMPLETE KIT** for ease, speed and cost effective use of resources.

Description	Size	Order Code
<b>Rocket HSG Procedure Pack</b> HSG Catheter Length 280mm, Maximum balloon volume: 3ml. Sterile, for single use, tin cartons of 10 procedure pack trays.	5FG	<b>R57422-00-PK</b>
	7FG	<b>R57424-00-PK</b>
<b>Rocket HSG Procedure Pack with integral wire stylet.</b> HSG Catheter length 280mm, Maximum balloon volume: 3ml. Sterile, for single use, in cartons of 10 procedure pack trays	5FG	<b>R57422-ST-PK</b>
	7FG	<b>R57424-ST-PK</b>