



Rocket® IPC offers an alternative long term method of peritoneal ascites drainage, allowing patients to be managed in the community which offers significant advantages:

- Early discharge
- Minimised repeat admission.
- Improved quality of life
- Improved control over discomfort and breathlessness
- Reduced cost

The system features:

PASSIVE DRAINAGE SYSTEM:

- Comfortable fluid flow, minimising discomfort during drainage.

SILICONE CATHETER:

- Soft, kink resistant material, giving increased patient comfort.
- 16Fg x 40cm catheter with 24cm fenestrated section, large diameter holes to maximise drainage and reduce occlusion risk.
- Complete with polyester cuff to promote rapid tissue in-growth, minimising infection risk and aiding catheter security.

SAFETY VALVE WITH UNIQUE LOCKING SYSTEM:

- Prevents leakage of air or fluid.
- Prevents accidental disconnection during drainage.
- Prevents accidental loss of valve caps, minimising contamination.

SIMPLE ORDERING:

- All the items needed for the drainage procedure and subsequent re-dressing are provided in a single pack

Description	Product Code
Rocket® IPC Dressing pack & drainage bag 2000ml drainage bag, gloves (latex free), field sheet, valve cap, foam catheter pad, 8 x gauze pads, clear-film dressing, 4 x alcohol swabs, blue clip, 'Tiger Striped' waste bag. Sterile for single use, in cartons of 5 packs	R54401 CH2 1968449
Rocket® Rocket IPC Valve Cap Sterile, for single use in cartons 10 caps	R54410-00-CP CH2 1968463
Rocket® Rocket IPC Drainage Line Sterile for single use, in cartons 10 lines	R54410-00-DL CH2 1968470
Rocket® IPC Replacement Valve Kit Sterile for single use, supplied individually	R54410-00-VV CH2 1968487