



**Rocket® IPC** is indicated for

- palliation of dyspnoea due to pleural effusion.
- pleurodesis as a treatment for pleural effusion

**Rocket® IPC** offers an alternative long term method of managing pleural drainage, allowing patients to be managed in the community which offers significant advantages:

- Early discharge
- Minimised repeat admission.
- Improved quality of life
- Improved control over discomfort and breathlessness
- Reduced cost

The system features:

**UNIQUE FLOW CONTROL:**

- Gives patient or carer total control over fluid flow, significantly reducing discomfort during drainage.
- Reduces risk of accidental vacuum release prior to catheter connection.

**SAFETY VALVE WITH UNIQUE LOCKING SYSTEM:**

- Prevents leakage of air or fluid.
- Prevents accidental disconnection during drainage.
- Prevents accidental loss of valve caps, minimising contamination.

Description	Product Code
<b>Rocket® IPC 600ml Dressing Pack &amp; Bottle</b> 600ml drainage bottle, gloves (non-latex), field sheet, valve cap, foam catheter pad, Gauze pads x 8, clear-film dressing, alcohol swabs x 4, blue closure clip, 'Tiger striped' waste bag. Sterile, for single use, in cartons of 5 packs	<b>R51400</b>
<b>Rocket® IPC 1000ml Dressing Pack &amp; Bottle</b> 1000ml drainage bottle with empty port, gloves (non-latex), field sheet, valve cap, foam catheter pad, Gauze pads x 8, clear-film dressing, alcohol swabs x 4, blue closure clip, 'Tiger striped' waste bag. Sterile, for single use, in cartons of 5 packs	<b>R51411</b>
<b>Rocket® IPC 600ml Mini Dressing Pack &amp; Bottle</b> 600ml drainage bottle, gloves (non-latex), field sheet, valve cap, alcohol swabs x 4. Sterile, for single use, in cartons of 5 packs	<b>R51400-00-PP</b>
<b>Rocket® IPC 1000ml Mini Dressing Pack &amp; Bottle</b> 1000ml drainage bottle with empty port, gloves (non-latex), field sheet, valve cap, alcohol swabs x 4. Sterile, for single use, in cartons of 5 packs	<b>R51411-00-PP</b>
<b>Rocket® IPC 1000ml Pre Evacuated Bottle</b> 1000ml drainage bottle, valve cap. Sterile, for single use, in cartons of 5 packs	<b>R51411-00-00</b>
<b>Rocket® Rocket IPC Valve Cap</b> Sterile, for single use in cartons 10 caps	<b>R54410-00-CP</b>
<b>Rocket® Rocket IPC Drainage Line</b> Sterile for single use, in cartons 10 lines	<b>R54410-00-DL</b>
<b>Rocket® IPC Replacement Valve Kit</b> Sterile for single use, in individual cartons	<b>R54410-00-VV</b>

