




- **KIT INCLUDES NRFit™ SPINAL NEEDLE:** 22G x 9cm QUINKE spinal needle
- **REDUCES RISK:** Combining all the necessary NRFit™ compatible components into a single kit significantly reduces the risk of mismatch or accidental connection to non-approved systems.
- **CLEAR, STURDY CLEAR CONSTRUCTION** for reliable, secure assembly and ease of use. Ready for immediate use for when time is precious.
- **TWO PART MANOMETER ASSEMBLY** allows CSF pressure monitoring between 3 - 40cmH₂O. Supplied with 3 way stopcock attached for connection to the NRFit™ compliant spinal needle and syringe included in the kit.



NRFit™

CONTENTS:

- 1 x NRFit™ compatible 40cm stackable Spinal Manometer Set
- 1 x NRFit™ compatible 22G x 9cm QUINKE spinal needle
- 1 x 10ml NRFit™ compatible syringe
- 5 x 10cm x 10cm gauze swabs
- 1 x 60ml gallipot
- 1 x 8" skin prep sponge stick
- 1 x fenestrated drape
- 1 x square drape 1m²
- 3 x 30ml Universal Containers with labels
- 1 x 20ml Syringe luer lock
- 1 x Microlance 21G /1.5" green needle (luer)
- 1 x Microlance 25G/5/8" yellow needle (luer)

Ordering Information
Code

Rocket® NRFit™ Lumbar Puncture Procedure Kit
complete ward or department procedure kit for single patient use,
packed in cartons 5 units

R55995-NRFIT

NRFit™ is a trademark of GEDSA and used with permission.