

with **NRFit™** connectors

ROCKET® SPINAL MANOMETERS with **NRFit™** connectors throughout the intrathecal fluid path to ensure connection to ISO80369-6 compliant spinal needles and syringes.


STACKABLE ASSEMBLY: **NRFit™** Male/female fittings of each tube permit additional sections to be stacked to allow pressures >40cmH₂O to be measured when required.

CLEAR, STURDY PLASTIC CONSTRUCTION for reliable, secure assembly and ease of use.

SUPPLIED COMPLETE with 3 way stopcock for attachment to all **NRFit™** ISO80369-6 compliant spinal needles and syringes.

READY FOR IMMEDIATE USE for when time is precious.

Description	Order Code	NHSSC Code
Rocket® Spinal Manometer Set with NRFit™ connectors Sterile, for single use, packed in cartons 10 units	R55990-ISO-6	eDirect GBR000



Rocket® Spinal Manometers have been providing safe, reliable, cost effective CSF pressure measurement for over 30 years.

The unique 'stackable' design, original to Rocket, has become a standard neurological diagnostic tool and can be found in Emergency Departments, Neurological, Medical Wards and specialist Neurosurgery centres throughout the world.

This original Rocket® design is shortly to be available with ISO80369-6 compliant connectors.

NRFit™ is a trademark of GEDSA and used with permission.