

## with NRFit™ connectors

**ROCKET® SPINAL MANOMETERS with NRFit™** connectors throughout the intrathecal fluid path to ensure connection to ISO80369-6 compliant spinal needles and syringes.


**STACKABLE ASSEMBLY:** NRFit™ Male/female fittings of each tube permit additional sections to be stacked to allow pressures >40cmH<sub>2</sub>O to be measured when required.

**CLEAR, STURDY PLASTIC CONSTRUCTION** for reliable, secure assembly and ease of use.

**SUPPLIED COMPLETE** with 3 way stopcock for attachment to all NRFit™ ISO80369-6 compliant spinal needles and syringes.

**READY FOR IMMEDIATE USE** for when time is precious.

Description	Order Code	NHSSC Code
<b>Rocket® Spinal Manometer Set with NRFit™ connectors</b> Sterile, for single use, packed in cartons 10 units	<b>R55990-ISO-6</b>	<b>eDirect GBR000</b>



Rocket® Spinal Manometers have been providing safe, reliable, cost effective CSF pressure measurement for over 30 years.

The unique 'stackable' design, original to Rocket, has become a standard neurological diagnostic tool and can be found in Emergency Departments, Neurological, Medical Wards and specialist Neurosurgery centres throughout the world.

This original Rocket® design is now available with ISO80369-6 NRFit™ compliant connectors.

NRFit™ is a trademark of GEDSA and used with permission.