



- **OVERCOME THE PROBLEMS** of high cost, cross infection risks and health and safety hazards commonly associated with reprocessing reusable LLETZ electrodes.
- **Rocket® LLETZ LOOPS** are designed to deliver higher cutting performance with less thermal damage to the sample.
- **2.4mm (3/32") x 10cm shaft** to fit all standard UK, US and European diathermy systems
- **DELTA™** style loops available for increased wire rigidity and increased control of penetration.

Rocket® LLETZ Range	Size		USA
<b>Rocket® LLETZ Diathermy Loop,</b> for single use, packed sterile in cartons 25 units	Small	14mm	<b>857800</b>
	Medium	18mm	<b>857801</b>
	Large	22mm	<b>857802</b>
	Square	10mm	<b>857803</b>
<b>Rocket® DELTA™ LLETZ Diathermy Loop,</b> for single use, packed sterile in cartons 25 units	Medium	18mm	<b>857851</b>
	Large	22mm	<b>857852</b>
<b>Rocket® LLETZ Ball Electrodes,</b> for single use, packed sterile in cartons 25 units	Standard	5mm	<b>857804</b>



North Island: 0800 658 814    South Island: 0800 804 546  
 customerservices@uslmedical.co.nz

**Rocket® LLETZ Loops** were the world's first single use loop electrodes designed specially for LLETZ.

**Rocket® LLETZ Loops** have been in continuous production for more than 27 years.

They continue our tradition of providing exceptional value, reliability and quality