

**ROCKET SPINAL MANOMETERS** with **LUER** connectors throughout the intrathecal fluid path to ensure connection to standard luer-lock spinal needles and syringes.

**READY FOR IMMEDIATE USE** for when time is precious.

**CLEAR, STURDY PLASTIC CONSTRUCTION** for reliable, secure assembly and ease of use.

**STACKABLE ASSEMBLY:** Male/female luer fittings of each tube permit additional sections to be stacked to allow pressures  $>40\text{cmH}_2\text{O}$  to be measured when required.

**SUPPLIED COMPLETE** with 3 way stopcock for attachment to all luer fitting spinal needles and syringes.

**Rocket® Spinal Manometers** have been providing safe, reliable, cost effective CSF pressure measurement for over 30 years.

The unique 'stackable' design, original to Rocket, has become a standard neurological diagnostic tool and can be found in Emergency Departments, Neurological, Medical Wards and specialist Neurosurgery centres throughout the world.



Description	Order Code	NHS SC Code
<b>Rocket Spinal Manometer Set</b> 1 x Rocket Spinal Manometer Set with 3 way stopcock.  Sterile, for single use, packed in cartons 40 units	<b>R55990</b>	<b>FTP002</b>