

ROCKET SPINAL MANOMETERS with **LUER** connectors throughout the intrathecal fluid path to ensure connection to standard luer-lock spinal needles and syringes.

READY FOR IMMEDIATE USE for when time is precious.

CLEAR, STURDY PLASTIC CONSTRUCTION for reliable, secure assembly and ease of use.

STACKABLE ASSEMBLY: Male/female luer fittings of each tube permit additional sections to be stacked to allow pressures $>40\text{cmH}_2\text{O}$ to be measured when required.

SUPPLIED COMPLETE with 3 way stopcock for attachment to all luer fitting spinal needles and syringes.

Rocket® Spinal Manometers have been providing safe, reliable, cost effective CSF pressure measurement for over 30 years.

The unique 'stackable' design, original to Rocket, has become a standard neurological diagnostic tool and can be found in Emergency Departments, Neurological, Medical Wards and specialist Neurosurgery centres throughout the world.



Description	Order Code	NHS SC Code
Rocket Spinal Manometer Set 1 x Rocket Spinal Manometer Set with 3 way stopcock. Sterile, for single use, packed in cartons 40 units	R55990	FTP002