



# Rocket Multipurpose Drainage Catheters

- **EasyTip™ CORE** to prevent 'ruckle-back', minimising resistance during insertion; reducing discomfort and speeding correct placement (8.4Fg-14FG).
- **SOFT, KINK FREE** catheter material, improves patient comfort and prevents unnecessary catheter replacement caused by harder, 'kink resistant' materials.
- **COMPLETE WITH BAG CONNECTION SET** to reduce costs and make Rocket® MPD Catheters even better value for your department or ward.
- **SECURE LOCKING SYSTEM** to give reliable pigtail locking without cumbersome 'keys' or bulky hubs. Self-adhesive removal instructions (for addition to the patient's notes) are provided with every locking drain.



**Rocket MPD Multipurpose Drainage Catheters and Sets offer a new solution to the challenge of maintaining a soft, kink free material with the need to prevent 'ruckle-back' during insertion.**



Description	Size	Code
<b>Rocket MPD Set: plain,</b> pigtail catheter, trocar and stiffener set, 30cm x 6 drainage holes complete with bag connection set	6FG	R58810-06-30
	8FG	R58810-84-25
	8FG	R58810-84-30
	14FG	R58810-14-30
<b>Rocket MPD Set: plain coated,</b> pigtail catheter, trocar and stiffener set, 30cm x 6 drainage holes complete with bag connection set	8FG	R58820-84-30
<b>Rocket MPD Set: locking pigtail catheter,</b> trocar and stiffener set, 30cm x 6 drainage holes complete with bag connection set	8FG	R58830-84-30
	10FG	R58830-10-30
<b>Rocket MPD Exchange Pigtail Catheter</b> 30cm x 6 drainage holes, plain	8FG	R58851-84-30