

Rocket Pro-Lux™ Proctoscope



The Rocket Pro-Lux™ Proctoscope is a fully disposable, ready to use, self-illuminated endoscopy system designed with new features to make usage convenient, practical and highly cost effective.

READY TO USE: self-contained battery powered light source with high intensity white light LED for maximum illumination.

NO EXTERNAL LIGHT SOURCE required for speed, simplicity and cost effective diagnosis and treatment.

COMPLETE SET supplied with a disposable obturator.

20mm INTERNAL DIAMETER for clear, precise vision and easy application of instrumentation. Ideal for Colorectal Out Patient or GP clinic use.



SAFETY: Should local hospital policy require it, the battery pack is fully removable to permit safe, environmentally responsible disposal. Each carton is supplied with a battery storage tray to ensure safe separation of used battery cells in conformance with COUNCIL DIRECTIVE 2006/66/EC (as amended) with regards to EU battery disposal regulations



CLIFFORDHALLAMHEALTHCARE

1300 720 274 www.ch2.net.au

Ordering Information	Rocket Code	CH2 Code
Rocket Pro-Lux™ Proctoscope Self-contained illumination, Non-sterile, 50 units per carton.	R53905	2129941

