




- **READY TO USE:** self-contained battery powered light source with a ring of six high intensity white light LED's for maximum illumination.
- **MAGNIFYING LENS** for clear and precise vision and easy application of biopsy instrumentation
- **COMPLETE SET** supplied with a disposable insufflation bulb and pilot balloon to provide a complete and economical single use system.
- **SAFETY:** Should local hospital policy require it, the battery pack is fully removable to permit safe, environmentally responsible disposal. Each carton is supplied with a battery storage tray to ensure safe separation of used battery cells in conformance with COUNCIL DIRECTIVE 2006/66/EC (as amended) with regards to EU battery disposal regulations.

The Rocket® Pro-Lux™ fully disposable sigmoidoscopy system is a ready to use device equipped with new features to make usage convenient, practical and effective.

The system utilises the well respected Welch-Allyn rigid sigmoidoscope combined with several additional features.



Description	NHS	Order Code
Rocket® Pro-Lux™ Sigmoidoscope Self-contained illumination handle, insufflation bulb, pilot balloon. Non-sterile, 25 units per carton	FFE1081	R53901
Rocket® Pro-Lux™ Sigmoidoscope - Lubricated Self-contained illumination handle, insufflation bulb, pilot balloon. Non-sterile, 25 units per carton.	FFE1078	R53901-L

nusurgix is a trading name of Rocket Medical plc