


- ◆ **LUER LOCK NEEDLE** provides up to a maximum of 10mm tissue penetration for increased safety.
- ◆ **SAFETY CLIP** prevents accidental deployment of the needle tip



- ◆ **REMOVABLE FINGER RING** can be re-sited at Position 1 or left as supplied in Position 2 to give best individual control and comfort during use.



**Rocket® Pudendal and Paracervical Block Needles** provide a safe and controllable method for delivering a pudendal or paracervical nerve block.

This route typically provides satisfactory perineal anaesthesia for normal delivery, low forceps manipulation and episiotomy.

The advantage of this approach over a major conduction block is the lack of circulatory and respiratory disturbance.

Typically the course and duration of labour are unaffected, even though there may be some impairment of the perineal-diaphragmatic-abdominal reflex.

Where there is concern about fetal compromise, either tested effective anaesthesia or, if time does not allow this, a pudendal block combined with local anaesthetic to the perineum can be used during instrumental birth.

Ref: Intrapartum care for healthy women and babies. NICE Clinical guideline. CG190 December 2014

Description	Order Code
<b>Rocket® Pudendal &amp; Paracervical Block Needles</b> 20G x 130mm. Sterile, for single use, individually packed, 20 per carton.	<b>R57040</b> <b>USL: 857040</b>



North Island: 0800 658 814 South Island: 0800 804 546  
[customerservices@uslmedical.co.nz](mailto:customerservices@uslmedical.co.nz)