



- ◆ **CE MARKED, PURPOSE DESIGNED DEVICE**, designed specifically for paracentesis and drainage of abdominal ascites, urinary bladder and other drainage applications which avoids the risks of inappropriate use of non-specific devices such as generic suprapubic drains.
- ◆ **SOFT, KINK FREE 8FG** catheter material, improves patient comfort and prevents unnecessary catheter blockage caused by thin walled materials.
- ◆ **SOFT INTEGRAL STITCH PLATE** combined with flexible tube support, provides simple, secure fixation after insertion.
- ◆ **EasyTip™ CORE** combined with a long taper tip helps prevent 'ruckle-back' and reduce resistance during insertion, so minimising discomfort and easing correct placement.
- ◆ **SAFETY SCALPEL** included as standard to aid compliance with COUNCIL DIRECTIVE 2010/32/EU
- ◆ **REMOVABLE NEEDLE CAP** permits aspiration to aid correct placement during insertion.
- ◆ **DRAINAGE BAG CONNECTION SET** included to reduce additional costs and make insertion and connection to a wide range of drainage options quick and secure.

Rocket® Multipurpose Needle Drains & Ward Procedure Packs offer a purpose designed, cost-effective solution for routine fluid drainage procedures in department or ward settings.

Indications:

For the drainage of body fluids including ascites, urine and pleural fluid.

For insertion into, but not restricted to, the abdominal cavity, chest and urinary bladder using a 'direct stab' technique combined with aspiration and/or ultrasound guidance.

Description	Order Codes
Rocket Multipurpose Needle Drain Set 8FG x 20cm x 6 holes with integral stitch plate and insertion needle set. 30cm drainage bag connection set with 2 way stopcock, No.11 retractable safety scalpel. Sterile, for single use, Carton: 5 sets.	R58800-08-ND CH2: 1969118



CLIFFORDHALLAMHEALTHCARE

1300 720 274 www.ch2.net.au

Rocket Medical Pty.

Suite 209, 1 Katherine Street, Chatswood. NSW 2067. Australia
 T: +61 (0) 431 146448 E: aus@rocketmedical.com W: aus.rocketmedical.com

Rev: 11
 Date: 2016-04-13
 Copyright© 2010-16