




- **SX™ SPIRAL ECHO BEVEL DESIGN** with high resolution echogenic surface preparation to give smoother, easier tissue penetration with extreme ultrasound definition.
- **SHAPED FINGER GRIP** with bevel marker and enhanced surfacing to significantly increase comfort and control
- **TWO NEEDLE STYLES:**
 - Rocket® Embryon®: dual angle bevel with 15mm brushed echogenic band
 - Rocket® SX™: back sharpened bevel, 40mm high definition echo spiral
- **TWO TUBING SET STYLES:**
 - Standard: two limb tube set with luer fittings into clear silicone bung
 - Oxford: single limb tube set into clear silicone bung with 45° metal cannula connection (right)



Rocket® Spiral Echo Needles are the latest generation of high performance, high ultrasound definition oocyte aspiration sets from Rocket Medical plc.



Description	Dimensions	Tube Set	Order Code
Rocket® Embryon® Single Lumen Needles Sterile, in cartons of 10 sets	17G x 32cm x 70cm	Oxford	R57602
	17G x 32cm x 70cm		R57603
	17G x 32cm x 70cm		R57605-00-90
Rocket® SX™ Single Lumen Aspiration Needles Sterile, in cartons of 10 sets	17G x 33cm x 70cm	Standard	R57600-SX-70
	16G x 33cm x 90cm		R57601-SX-90
	17G x 33cm x 55cm	Oxford	R57602-SX-55
	17G x 33cm x 70cm		R57602-00-SX
	17G x 33cm x 90cm		R57602-SX-90
	17G x 35cm x 70cm		R57603-SX-70
	17G x 35cm x 90cm		R57603-SX-90