


 URIPOT[®]

 URIKONE[®]

Both designs feature:



- **UNIVERSAL CONTAINER**, used as standard by most laboratory services, makes urine collection easy and cost effective.
- **PRE-ASSEMBLED**: simple and ready to use, suitable for either male or female.
- **INCLUDES AQUEOUS PRE WIPE** for easy, rapid cleansing of the genital area prior to sampling.
- **PAEDIATRIC & NEONATAL USE**: Both Uripot[®] and Urikone[®] are suitable for paediatric and neonatal collection when used under nursing or parental supervision to ensure adequate sample collection and child safety.

Rocket[®] URIKONE[®] & URIPOT[®] offer cost effective, ready to use urine collection in two styles to meet the widest range of patient needs.

URIPOT[®] has an anatomically designed filling pot with soft, smooth edges and is warm to the touch, making sampling easier and more comfortable for the female patient.

A special non return valve makes gradual filling simple and helps prevent overflow and spillage.

URIKONE[®] offers a simple, traditional cone design ideal for male urine collection.

Description		Code
 Supplied individually wrapped, sterile, in cartons of 50 units	with prep wipe	R53084
	without wipe	R53086
	with boric acid	R53087
	with boric acid & wipe	R53088
 Supplied individually wrapped, sterile, in cartons of 50 units	with prep wipe	R53100
	without wipe	R53102