

Rocket Water Trap Set for use with Rocket Craft™ Suction Pumps

The Rocket Water Trap Set for use with all Rocket Craft™ Suction Pumps offers a simple, cost effective solution to contamination risks and reuse problems.

- **COMPLETE REPLACEMENT SET:** Provides the complete bottle and tubing set in single, sterile pack for one time patient use.
- **GLASS TRAP BOTTLE:** To maintain clear view in the event of potential pump contamination with aspirate. Designed to withstand >450mmHg vacuum.
- **REMOVES CONTAMINATION RISK** from irregular or uncontrolled decontamination of reusable bottle, bung and angle tube set.
- **PATIENT TUBE SET:** 2m patient tube set with integral luer lock connector to ensure secure, reliable connection to the aspiration needle set.
- **SUITABLE FOR USE:** Compatible with all versions of Rocket® Craft™ Suction and Duo Vac. Also compatible with the *genX and Pioneer Pro Pumps*.



Ordering Information

Rocket Water Trap Set:

suitable for all Craft™ Suction Pumps packed individually, for single patient use, 10 per carton

R57685



Rocket Medical plc 150 Recreation Park Drive Unit 1 Hingham MA 02043 USA
T: 1-781-749-6223 **F:** 1-781-749-6235 **E:** usa@rocketmedical.com www.rocketmedical.com