

Rocket

Digital Oocyte Aspiration Pump Filter Sets

Rocket Digital Oocyte Aspiration Pump Filter Sets offer a simple, cost effective, patient connection set specifically designed for use with the latest Rocket Digital Oocyte Aspiration Pump.

- **REMOVES CONTAMINATION RISK** from irregular or uncontrolled decontamination of multiple use consumable filters.
- **PATIENT TUBE SET:** 2.5m lightweight patient tube set with integral luer lock connector to ensure secure, reliable connection to all standard luer-lock fitting aspiration needle sets.
- **HYDROPHOBIC BARRIER** provides pump protection by preventing accidental aspiration of fluid into the device.
- **DESIGNED SPECIFICALLY** to match the performance characteristics of the Rocket Digital Oocyte Aspiration Pump



FILTER SPECIFICATION

17mm
 Effective filtration 0.2micron,
 BFE 99.9999%*
 VFE 99.9999%*
 High flow /low resistance Insufflation/Gas/Air Filter
 Hydrophobic barrier.
 Inner core: PTFE

Description	Code
Rocket Digital Oocyte Aspiration Pump (Patient Connection) Filter Set 0.2µm filter set, 2.5m patient connection tube with male luer connector. Supplied sterile, for single patient use in cartons 10 units.	R57686
Rocket Aspiration Pump Filter Set for COOK KMAR Series Pumps 0.2µm filter set, 2.5m patient connection tube with male luer connector. Supplied sterile, for single patient use in cartons 10 units.	R57687





Rocket Medical plc 150 Recreation Park Drive Unit 1 Hingham MA 02043 USA
T: 1-781-749-6223 **F:** 1-781-749-6235 **E:** usa@rocketmedical.com www.rocketmedical.com

Registered office: Rocket Medical plc. Imperial Way WATFORD WD24 4XX England. Company Number: 03276608