

ROCKET SPINAL MANOMETERS with LUER connectors throughout the intrathecal fluid path to ensure connection to standard luer-lock spinal needles and syringes.

READY FOR IMMEDIATE USE for when time is precious.

CLEAR, STURDY PLASTIC CONSTRUCTION for reliable, secure assembly and ease of use.

STACKABLE ASSEMBLY: Male/female luer fittings of each tube permit additional sections to be stacked to allow pressures $>40\text{cmH}_2\text{O}$ to be measured when required.

SUPPLIED COMPLETE with 3 way stopcock for attachment to all luer fitting spinal needles and syringes.

Rocket® Spinal Manometers have been providing safe, reliable, cost effective CSF pressure measurement for over 30 years.

The unique 'stackable' design, original to Rocket, has become a standard neurological diagnostic tool and can be found in Emergency Departments, Neurological, Medical Wards and specialist Neurosurgery centres throughout the world.



| Description | Order Code |
|---|---------------|
| Rocket Spinal Manometer Set 1 x Rocket Spinal Manometer Set with 3 way stopcock. Sterile, for single use, packed in cartons 40 units | R55990 |