



- **SURGICAL GRADE STAINLESS NEEDLE** and 'dual eyed' sprung obturator are securely moulded into a specially shaped hub.
- **DUAL OPPOSED OBTURATOR 'EYES'** reduces risk of blockage from tissue and ensures rapid fill so reducing operating time.
- **INJECTION MOULDED HUB** assembly ensures strength and needle security with no risk of needle detachment.
- **COMFORTABLE, EASY GRIP HAND PIECE** gives precise control on insertion
- **UNIQUE RED/GREEN INDICATOR SYSTEM** shows exact position of the obturator ensuring the needle point is not unknowingly exposed
- **INSUFFLATION FILTER SETS** with 0.1 micron filtration and hydrophobic barrier complete the range.



**Rocket® Verres™ Insufflation Needles** offer a safe and cost effective method of gas insufflation for laparoscopy avoiding difficult and expensive reprocessing which is so frequently the cause of blunt, bent or blocked needles and jammed spring mechanisms.

Description	Size/Length	Code
<b>Rocket® Verres™ Insufflation Needle</b> Carton 20 units, sterile for single patient use	14g x 120mm	<b>R57920</b> <b>USL: 857920</b>
	14g x 150mm	<b>R57921</b> <b>USL: 857921</b>



North Island: 0800 658 814    South Island: 0800 804 546  
customerservices@uslmedical.co.nz