


KleenSpec[®]

Disposable Sigmoidoscopes & Accessories



Rocket Medical are dedicated to the design and manufacturer of their own products and very rarely would we would consider distribution for another producer – Welch-Allyn UK is one such company. They have a long and distinguished tradition of quality and innovation and we are delighted to be their distribution partner in the UK.



- 19 mm x 25 cm supplied complete with obturator, Available plain or lubricated, with or without filters
- High-visibility graduations and glare free frosted coating for maximum light output at tip
- Hydrophobic filters prevent internal contamination of the inflation bulb
- Disposable for economy, safety and convenience
- Supplied non-sterile, individually wrapped in carton 100 units.

Welch Allyn Code	Rocket order code	Welch Allyn Code	Rocket order code
53130 Disposable Sigmoidoscope ① unlubricated,	R53130	36019 Light Head for disposable ② sigmoidoscopes	R26360-19
53130-L Disposable Sigmoidoscopes ① lubricated	R53130-L	73211 Light handle with lead ④	R26732-11
53130-F Disposable Sigmoidoscopes ① unlubricated with filters	R53130-F	73324 Transformer 6V with lead ⑥	R26733
53130-LF Disposable Sigmoidoscopes lubricated with filters	R53130-LF	30200 Insufflation Bulb Set, ⑤ Complete	R26302
30210 Pack of 10 Hydrophobic Filters ③	R53302-10	73500 3.5V Halogen Handle Adaptor	R26735
30215 Luer Lock for use with 30210	R26302-15		

KleenSpec[®] and the Welch Allyn[®] logo are registered trademarks of Welch Allyn UK, other logos and designs are trademarks of Rocket Medical plc.

