

Single Chest Drain Bottle

R54500 & R54509

Dear Customer

As part of our programme of continuous improvement we are pleased to announce an update for the R54500 & R54509 Chest drainage bottles.

Rocket Medical have developed a droplet control shroud.

The shroud when used in conjunction with the R54571 Filter and R54569 Portable suction unit will allow the PSU to be used on patients with a viral infection such as Covid-19.



Customer Information Bulletins on all of our products can be found on our web site at:
www.rocketmedical.com

If you have any questions, please contact the Customer Services Team for assistance:

Tel: +44 191 419 6988. Fax: +44 419 6889. Email: customerservices@rocketmedical.com