

# CUSTOMER INFORMATION BULLETIN

## Rocket Information for Use Updates

Dear Customer,

Please be aware that the Instructions for Use (IFU) have been updated to remove reference to the use of fluoroscopy upon insertion of the following devices:

Product Code (IFU Ref)	Product Description
<b>R58800-08-SD</b> (ZDOCK237)	Rocket Safety Needle Drain 8FG x 16cm paracentesis catheter with blunt obturator needle
<b>R58800-06-BD</b> (ZDOCK179)	Multipurpose Drain with fixation plate 6Fg x 20cm Includes insertion needle and bag connection set
<b>R58800-06-BK</b> (ZDOCK179)	Multipurpose Drain procedure Kit with fixation plate, 6fg x 20cm, Complete ward drainage procedure pack

This update has been made following the review of guidance<sup>1,2</sup> that indicates the use of Ultrasound during paracentesis is best practice.

**Please note this update to the IFU does not affect patient safety.**

#### References:

1. Aithal GP, et al. Guidelines on the management of ascites in cirrhosis. *Gut* 2020;0:1–21.
2. Recommendations on the Use of Ultrasound Guidance for Adult Abdominal Paracentesis: A Position Statement of the Society of Hospital Medicine. 2019; 14: E7–E15

Thank you,



Richard Varey

Research & Development Product Manager

Rocket Medical

Customer Information Bulletins on all of our products can be found on our web site at:

[http://www.rocketmedical.com/Rocket\\_support\\_CIB\\_main.htm](http://www.rocketmedical.com/Rocket_support_CIB_main.htm)

**If you have any questions, please contact the Customer Services Team for assistance:**

Tel: +44 (0) 191 419 6988. Fax: +44 419 6889. Email: [customerservices@rocketmedical.com](mailto:customerservices@rocketmedical.com)



Rocket Medical plc, Sedling Road Washington Tyne & Wear NE38 9BZ England

T: +44 (0) 191 419 6988 F: +44 (0) 191 419 6989 E: [customerservices@rocketmedical.com](mailto:customerservices@rocketmedical.com) [www.rocketmedical.com](http://www.rocketmedical.com)

Registered Office: Rocket Medical plc, Imperial Way, Watford, WD24 4XX, England Company No: 03276608