

Rocket® EZI Plug Fetal Scalp Electrode Discontinuation

Dear Customer,

We recently communicated the obsolescence of our EZI plug Fetal Scalp Electrode product range. In order to give you greater time to transition to an alternative, we will now cease supply of the Rocket® EZI plug FSE on 31st August, 2020. The following product code is affected:

Stock code	NHS Code	Description
R57003-EZI	FDK101	Rocket® EZI plug FSE

The recommended alternative is the Rocket® Copeland FSE. This has a similar design to the Rocket® EZI FSE and is available with a Quik connector for Corometrics, Huntleigh and Hewlett Packard systems or with a connector for the Philips Avalon FM30 and FM50 systems.

Stock code	NHS Code	Description
R57008-00-CN	FDK097	Rocket Copeland Fetal Scalp Electrode with Quik connector for Corometrics, Huntleigh and Hewlett Packard
R57008-00-FM	FDK557	Rocket Copeland Fetal Scalp Electrode with Kwik connector for Philips Avalon FM30 and FM50

Please note that in order to use these electrodes you will also need to change your leg plate cable and leg plate electrodes. The relevant ordering codes for these items for Huntleigh monitors are shown below. To order these items or further information please contact Huntleigh.

For information about compatible cables and leg plate electrodes for other monitors, please contact your monitor manufacturer.

Huntleigh Monitoring System	Leg Plate Cable Required	Leg Plate Electrode Required
Team 3	ACC-OBS-068	ACC-2464AA0
FM800 Encore	ACC-OBS-023	ACC-2464AA0
Team (IP) and FM800 Adapt	ACC-739087	ACC-2464AA0

If you have recently purchased the Rocket® EZI fetal scalp electrodes, a Rocket Medical representative will be in touch soon to discuss how we can support you during this transition. If, in the meantime, you have any questions please contact the Customer Services Team for assistance:

Tel: +44 (0) 191 419 6988

Email: customerservices@rocketmedical.com

Customer Information Bulletins for all of our products can be found on our website at: www.rocketmedical.com