

# CUSTOMER INFORMATION BULLETIN

## Rocket® R54500 & R54509 Single Chest Drainage Bottle

Dear Customer,

As part of our programme of continuous improvement we are pleased to introduce an update for the R54500 & R54509 Chest Drainage bottles.

- There is now an automatic positive pressure relief valve.
- If the flow of escaping air is blocked, the positive pressure relief valve will allow the air to escape and protect the patient from potential harm.



**Positive Pressure Blow Off Valve**

Customer Information Bulletins on all of our products can be found on our web site at: [www.rocketmedical.com](http://www.rocketmedical.com)

**If you have any questions, please contact the Customer Services Team for assistance:**

Tel: +44 191 419 6988. Fax: +44 419 6889. Email: [customerservices@rocketmedical.com](mailto:customerservices@rocketmedical.com)

