

Dear Customer,

Following a review of product R55995-NRFIT® Lumbar Puncture Tray with NRFit® manometer, we have taken the decision to remove the Pajunk QUINKE NRFit® 22G x 9cm spinal needle. Customers have told us that they have an individual preference for the type of spinal needle used; removing the needle allows users to select a NRFit® needle of their choice.

REVISED KIT CONTENTS:

- 1 x NRFit® 40cm H₂O stackable spinal manometer
- 1 x 10ml NRFit® slip syringe
- 5 x 10cm x 10cm gauze swabs
- 1 x 60ml pot
- 1 x Skin prep sponge stick
- 1 x 2-piece fenestrated drape,
75 x 90cm (75 x 45cm x 2)
- 1 x square drape 60 x 60 cm
- 3 x 30ml Universal Containers with labels
- 1 x 20ml luer lock syringe
- 1 x Microlance 21G /1.5" green
needle (luer)
- 1 x Surguard 2 25G /38mm orange
needle (luer)



Ordering Information	Code
Lumbar Puncture Tray with NRFit® Spinal Manometer Complete procedure kit, for single patient use. Packed in cartons of 5 units.	R55995-NRFIT

If you have any questions, please contact the Customer Services Team or your local Rocket Medical representative for assistance:

Tel: +44 (0) 191 419 6988

Email: customerservices@rocketmedical.com

Customer Information Bulletins for all of our products can be found on our website at:

www.rocketmedical.com

NRFit® is a trademark of GEDSA and used with permission.