



Dear Customer

We have recently been approached for advice regarding the use of the Filter material used in the Rocket Medical PSU Filter and the Rocket Medical Pleural Vent on COVID-19 patients and the potential for the control of aerosols.

The filter materials used are designed for the control of fluid and are not viral filters.

Following advice from a number of leading clinicians, we understand that hospitals are treating COVID-19 patients in dedicated areas. In this case the Pleural Vent or under water seal chest drains could be used as normal with appropriate PPE.

Additionally, if a greater level of assurance is required, connecting the patient to thoracic wall suction with an appropriate filter may reduce any perceived risk of infection.

Customer Information Bulletins for all of our products can be found on our web site at:
www.rocketmedical.com

If you have any questions, please contact the Customer Services Team for assistance:
Tel: +44 191 419 6988. Email: customerservices@rocketmedical.com