



Dear Customer

We have recently been approached for advice regarding the use of Under Water Seal (UWS) chest drains on COVID-19 patients and the potential for the release of aerosols.

Following advice from a number of leading clinicians, we understand that hospitals are treating COVID-19 patients in dedicated areas. In this case UWS chest drains could be used as normal with appropriate PPE.

Additionally, if a greater level of assurance is required, connecting the patient to thoracic wall suction with an appropriate filter may reduce any perceived risk of infection.

Customer Information Bulletins for all of our products can be found on our web site at:
www.rocketmedical.com

If you have any questions, please contact the Customer Services Team for assistance:
Tel: +44 191 419 6988. Email: customerservices@rocketmedical.com