



The system features:

- ◆ **UNDERWATER SEAL CHEST DRAINAGE BOTTLE** for either single collection (R54500) or double collection (R54509).
- ◆ **SAFETY PRESSURE RELIEF VALVE** provides a secondary escape route for air preventing positive pressurisation of the drainage bottle.
- ◆ **UPRIGHT, INHERENTLY STABLE DESIGN** which does not require tiresome stands or carriers.
- ◆ **1800ml CAPACITY:** available in two graduations 25ml (R54500) or 5ml (R54509).
- ◆ **AVAILABLE WITH A WIDE RANGE OF TUBING SETS:** including, single cone connection to large bore catheters: (R54502) and also for use with Rocket® Seldinger Chest Drainage Sets. (R54539).



**The Rocket® Blue Bottle™** single chamber chest drain and tube set combine to provide the UK's best selling chest drainage system.

The universally accepted technique of choice for pleural drainage following cardiac and thoracic surgery, in A&E departments and medical wards.

- Traditional underwater seal
- Simple, logical functionality
- Easy setup

**Function:** Traditionally, the simplest method to accomplish chest drainage, is to assemble a single 'bottle' containing a prescribed amount of water. One short tube leading to the outside atmosphere functions as a vent; whilst another leads from the patient and is submerged approximately 2cm below the surface of the water.

This creates an 'underwater seal', the most important element in pleural drainage. The water-seal provides a low resistance, one way valve that allows air and fluid to be expelled from the pleural space by means of positive expiratory pressure and prevents atmospheric air from being drawn back into the pleural cavity.

Description	Order Code
<b>Rocket® Blue Bottle™</b> Single chamber, 1800ml capacity Sterile for single use, 10 units per carton.	<b>R54500</b> <b>NHSFSW076</b>
<b>Rocket® Double Chamber bottle</b> ideal for paediatric or specific drainage recording. 1800ml capacity bottle, 1st chamber: 5ml, 2nd chamber: 10ml graduations, sterile, 10 units per carton	<b>R54509</b> <b>NHSFSW047</b>
<b>Rocket® Tubing Set</b> with fir tree connector with bottle cap, integral straw, support spring, 1.5mm , sterile, single wrapped, 10 units per carton	<b>R54502</b> <b>NHSFSW077</b>
<b>Rocket® Tubing Set</b> with plain ended connection suitable for connection to Rocket® Seldinger Chest Drain Sets, with bottle cap, integral straw, support spring, 1.5mm. Sterile for single use, 10 units per carton	<b>R54539</b> <b>NHSFSW187</b>