



- ◆ Safe, quick and easy to apply without the need for applicators.
- ◆ Ribbed finger grip for easier placement.
- ◆ Can be applied with as little as 1.5cm cervical dilation for earlier monitoring.
- ◆ Three point contact for excellent trace quality
- ◆ Automatic control of depth of penetration of needle for low scalp trauma.
- ◆ Cannot detach accidentally when correctly applied.
- ◆ Easy to reposition and remove.



The Rocket Copeland® name and Fetal Scalp Monitoring are synonymous and for over 20 years Copeland® Fetal Scalp Electrodes have provided safe and reliable recording of the fetal heart rate during labour.

Description	USL	Rocket
Rocket® Copeland® Fetal Scalp Electrode for Qwik Connect Cables Sterile, for single use, individually packed, 25 per carton.	857008-00-CN	R57008-00-CN
Rocket® Copeland® Fetal Scalp Electrode for Philips Avalon FM30/FM50 DECG Cables Sterile, for single use, individually packed, 25 per carton.	857008-00-FM	R57008-00-FM