

Rocket® Double Lumen Aspiration Needles Reproductive Medicine



- **FLUSHING ARM STOPCOCK** fitted as standard to minimise accidental aspiration of flushing media and increase flushing control.
- **SHAPED FINGER GRIP** with bevel marker and enhanced surfacing to significantly increase comfort and control
- **BEVEL STYLE:**
 - Rocket® Embryon® dual angle bevel with 15mm brushed echogenic band
- **OXFORD TUBING SET:** Single limb tube set into a clear silicone bung with a 45° metal cannula connection (as illustrated)

Description	Dimensions	Order Code
Rocket® Double Lumen Oocyte Aspiration Needles etched echo band, polyurethane, Oxford style tube set with clear silicone stopper for Falcon® tubes, pouch packed, supplied sterile, 10 per carton, for single patient use	16g x 32cm x 70cm	R57610
	16g x 35cm x 70cm	R57611
	16g x 35cm x 90cm	R57611-00-90