



- **OVERCOME THE PROBLEMS** of high cost, cross infection risks and health and safety hazards commonly associated with reprocessing reusable LLETZ electrodes.
- **Rocket[®] LLETZ LOOPS** are designed to deliver higher cutting performance with less thermal damage to the sample.
- **2.4mm (3/32") x 10cm shaft** to fit all standard UK, US and European diathermy systems.
- **DELTA™** style loops available for increased wire rigidity and increased control of penetration.

Rocket[®] LLETZ Loops were the world's first single use loop electrodes designed specially for LLETZ.

Rocket[®] LLETZ Loops have been in continuous production for more than 27 years. They continue our tradition of providing exceptional value, reliability and quality

Rocket [®] LLETZ Range	Size		US
Rocket[®] LLETZ Diathermy Loop, for single use, packed sterile in cartons 25 units	Small	14mm	857800
	Medium	18mm	857801
	Large	22mm	857802
	Square	10mm	857803
Rocket[®] LLETZ Ball Electrodes, for single use, packed sterile in cartons 25 units	Standard	5mm	857804
Rocket[®] DELTA™ LLETZ Diathermy Loop, for single use, packed sterile in cartons 25 units	Medium	18mm	857851
	Large	22mm	857852



North Island: 0800 658 814 South Island: 0800 804 546
customerservices@uslmedical.co.nz