



- ◆ **CE MARKED, PURPOSE DESIGNED DEVICE**, designed specifically for paracentesis and drainage of abdominal ascites, urinary bladder and other drainage applications which avoids the risks of inappropriate use of non-specific devices such as generic suprapubic drains.
- ◆ **SOFT, KINK FREE 8FG** catheter material, improves patient comfort and prevents unnecessary catheter blockage caused by thin walled materials.
- ◆ **SOFT INTEGRAL STITCH PLATE** combined with flexible tube support, provides simple, secure fixation after insertion.
- ◆ **EasyTip™ CORE** combined with a long taper tip helps prevent 'ruckle-back' and reduce resistance during insertion, so minimising discomfort and easing correct placement.
- ◆ **SAFETY SCALPEL** included as standard to aid compliance with COUNCIL DIRECTIVE 2010/32/EU
- ◆ **REMOVABLE NEEDLE CAP** permits aspiration to aid correct placement during insertion.
- ◆ **DRAINAGE BAG CONNECTION SET** included to reduce additional costs and make insertion and connection to a wide range of drainage options quick and secure.

**Rocket® Multipurpose Needle Drains & Ward Procedure Packs offer a purpose designed, cost-effective solution for routine fluid drainage procedures in department or ward settings.**

**Indications:**

For the drainage of body fluids including ascites, urine and pleural fluid.

For insertion into, but not restricted to, the abdominal cavity, chest and urinary bladder using a 'direct stab' technique combined with aspiration and/or ultrasound guidance.



**R58800-08-NK Rocket Multipurpose Ward Procedure Pack**

Description	Order Codes
<b>Rocket Multipurpose Needle Drain Set</b> 8FG x 20cm x 6 holes with integral stitch plate and insertion needle set. 30cm drainage bag connection set with 2 way stopcock, No.11 retractable safety scalpel. Sterile, for single use, Carton: 5 sets.	<b>R58800-08-ND</b>
<b>Rocket Multipurpose Needle Drain Procedure Pack</b> 8FG needle drain with trocar needle, 30cm drain bag connection set + 2 way tap, 2000ml drainage bag with empty port, 20ml, 12ml syringes, luer lock, 21g x 1.5" needle, 25g x 5/8" needle, 5in Sharp/Blunt Scissors, 60ml gallipot, Pair medium polycy gloves, gauze swabs x 5, No.11 retractable safety scalpel, 11x15cm 'Mepore' dressing, 50x60m adhesive fenestrated drape, Sterile, for single use, Carton: 5 trays	<b>R58800-08-NK</b>