


Band-It™ Haemorrhoid Ligator



The Band-It™ single use, multiband haemorrhoid ligator comes pre-loaded with five hypoallergenic ligature bands for single-handed deployment of up to five bands without the inconvenience of re-loading.

- **READY TO USE:** Self-contained, pre-loaded with five hypoallergenic haemorrhoid bands.
- **COMPLETE SET:** Supplied with soft, highly flexible suction tube for reduced drag.
- **INDEXED DEPLOYMENT WHEEL:** Makes band delivery quick and easy to control.
- **SIMPLE OPERATION:** Designed for single handed deployment with either the left or the right hand.

Description	Style	Order Code
Band-It™ Haemorrhoid Ligator Non-sterile, 20 units per carton	Hypoallergenic bands with suction tube	R53703 CH22211071


CLIFFORDHALLAMHEALTHCARE

1300 720 274 www.ch2.net.au