





Haemorrhoid Injection Sets provide a safe, reliable method of injection of Oily Phenol BP, or other sclerosing agents, as an outpatient treatment for haemorrhoids.



The design features:

- **COMFORTABLE** three finger ring grip for ease of injection of viscid Phenol.
- **ANGLED or STRAIGHT** luer lock needle for security and improved ease of access.
- **WIDE BORE** needle with 18G distal tip to provide maximum control during injection.
- **SUPPLIED COMPLETE** with filling cannula for ease of use in the busy clinic setting.
- **RELIABLE**, practical and competitively priced to ensure value for your department.

Ordering Information		
Haemorrhoid Injection Set Supplied with 10ml, 3 ring syringe, haemorrhoid injection needle, filling cannula. Supplied sterile, for single use. 25 sets per carton.	Straight, 10cm, 4"	856210
	Angled, 10cm, 4"	856211



North Island: 0800 658 814 South Island: 0800 804 546
 customerservices@uslmedical.co.nz