

# URIPOT™ Urine Collection Pack



- **UNIVERSAL CONTAINER**, used as standard by most laboratory services, makes urine collection easy and cost effective.
- **PRE-ASSEMBLED**: simple and ready to use, suitable for either male or female.
- **INCLUDES AQUEOUS PRE-WIPE** for easy, rapid cleansing of the genital area prior to sampling.
- **PAEDIATRIC & NEONATAL USE**: URIPOT™ is suitable for paediatric and neonatal collection when used under nursing or parental supervision to ensure adequate sample collection and child safety.

The Rocket Medical URIPOT™ offers a cost effective, ready to use urine collection in two styles to meet the widest range of patient needs.

URIPOT™ has an anatomically designed filling pot with soft, smooth edges and is warm to the touch, making sampling easier and more comfortable for the female patient.

A special non return valve makes gradual filling simple and helps prevent overflow and spillage.

Description	Description	Code
URIPOT™ Urine Collection Pack Supplied individually wrapped, sterile, in cartons of 50 units	with Wipe	<b>R53084</b> <b>USL: 853084</b>
	without Wipe	<b>R53086</b> <b>USL: 853086</b>



North Island: 0800 658 814    South Island: 0800 804 546  
customerservices@uslmedical.co.nz