


# Rocket<sup>®</sup> Amniotomy Hook



- ◆ **RE-DESIGNED** shaft and distal hook to improve rigidity and hook profile
- ◆ **GUARDED HOOK:** located at a controlled depth to permit perforation of the membranes without trauma to the surrounding tissue
- ◆ **CURVED:** anatomically shaped to ensure a positive grip in the users hand and to follow the normal curved pelvic outlet. for precise finger location and positioning.

Description	Order Code
<b>Rocket<sup>®</sup> Amniotomy Hook</b> , curved, 25cm. Sterile, for single use, individually packed, 100 per carton.	<b>R57065</b>