



Band-It™ Haemorrhoid Ligator



The Band-It™ single use multiband haemorrhoid ligator comes pre-loaded with five hypoallergenic ligature bands for single-handed deployment of up to five bands without the inconvenience of re-loading.

- **READY TO USE:** Self-contained, pre-loaded with five hypoallergenic haemorrhoid bands.
- **COMPLETE SET:** Supplied with soft, highly flexible suction tube for reduced drag.
- **INDEXED DEPLOYMENT WHEEL:** Makes band delivery quick and easy to control.
- **SIMPLE OPERATION:** Designed for single handed deployment with either the left or the right hand.

Description	Style	Rocket Code	
Band-It™ RBL Haemorrhoid Ligator Non-sterile, 20 units per carton	Hypoallergenic bands with suction tube	R53703	853703



North Island: 0800 658 814 South Island: 0800 804 546

customerservices@uslmedical.co.nz