


Design features:

- **UNIVERSAL CONTAINER**, used as standard by most laboratory services, makes urine collection easy and cost effective.
- **PRE-ASSEMBLED**: simple and ready to use, suitable for either male or female.
- **INCLUDES AQUEOUS PRE-WIPE** for easy, rapid cleansing of the genital area prior to sampling.
- **PAEDIATRIC & NEONATAL USE**: Both URIPOT™ and URIKONE™ are suitable for paediatric and neonatal collection when used under nursing or parental supervision to ensure adequate sample collection and child safety.
- **CONTAINS BORIC ACID** for microbiological analysis

URIPOT™ & URIKONE™ offer cost effective, ready to use urine collection in two styles to meet the widest range of patient needs.

URIPOT™ has an anatomically designed filling pot with soft, smooth edges and is warm to the touch, making sampling easier and more comfortable for the female patient.

A special non return valve makes gradual filling simple and helps prevent overflow and spillage.

URIKONE™ offers a simple, traditional cone design ideal for male urine collection.

Description		Code
URIPOT™ Urine Collection Pack Supplied individually wrapped, sterile, in cartons of 50 units	with Wipe	R53084
	without wipe	R53086
	with Boric Acid	R53087
	with Boric Acid & Wipe	R53088
URIKONE™ Urine Collection Pack Supplied individually wrapped, sterile, in cartons of 50 units	with Wipe	R53100
	without wipe	R53102