

HSG Procedure Packs



Rocket Medical HSG Procedure Packs provide the necessary single use items for cost effective uterine canalisation for hysterosalpingography.

CONTENTS:

- HSG Catheter
- Vaginal speculum
- Sponge stick for prep solution
- Sponge holder with 5 x green, x-ray detectable swabs
- Cervical dilator
- Water-soluble lubricant
- 4 x 30cm swab sticks
- Filling kwill
- 20ml and 3ml syringes
- 30cm contrast extension tube

Rocket Medical HSG Catheters have these well-established and trusted features:

- **LATEX-FREE** balloon minimises the risk of patient tissue reactions.
- **SMOOTH TIP AND LOW PROFILE BALLOON** ensures low trauma insertion and removal – ideal for the nulliparous cervix.
- **REMOVABLE STIFFENING SHEATH** to aid insertion in the heavily retroverted or antverted cervix.
- **CONTRAST PATH TAP** ensures retention of contrast when further volumes of media are required.
- **INTEGRAL WIRE STYLET** for improved malleability without using tiresome removable stylets (ST versions)
 - ◇ Improves ease of insertion
 - ◇ Increased value for money
- For use with aqueous contrast media.
- **COMPLETE KIT** for ease, speed and cost-effective use of resources.



CLIFFORDHALLAMHEALTHCARE

1300 720 274 www.ch2.net.au

Description	Size	Rocket Code	CH2 Code
Rocket HSG Procedure Pack Catheter length 280mm Maximum balloon volume: 2ml Sterile, for single use, in cartons of 10	5FG	R57422-00-PK	2176641
	7FG	R57424-00-PK	2176683
Rocket HSG Procedure Pack with Stylet Catheter length 280mm Maximum balloon volume: 2ml Sterile, for single use, in cartons of 10	5FG	R57422-ST-PK	2176691
	7FG	R57424-ST-PK	2176704