


HSG & HyCoSy Catheters



Rocket Medical HSG Catheters provide safe and effective uterine canalisation for hysterosalpingography and hysterosalpingo contrast sonography.

CONTENTS:

- HSG and HyCoSy Catheter
- Introducer sheath
- Pre-attached syringe
- Protective balloon cover

Rocket Medical HSG and HyCoSy Catheters have these well-established and trusted features:

- **LATEX-FREE** balloon anchors the catheter and provides a contrast seal.
- **SMOOTH TIP AND LOW PROFILE BALLOON** ensures low trauma insertion and removal – ideal for the nulliparous cervix.
- **REMOVABLE STIFFENING SHEATH** to aid insertion in the heavily retroverted or anteverted cervix.
- **CONTRAST PATH TAP** ensures retention of contrast when further volumes of media are required.
- **INTEGRAL WIRE STYLET** for improved malleability without using tiresome removable stylets (ST versions)
 - ◇ Improves ease of insertion
- **COMPATIBLE** with aqueous based contrast media.

| Description | Size | Order Code | CH2 Code |
|--|------|---------------------|----------------|
| HSG and HyCoSy Catheter HSG Catheter Length 280mm, Maximum balloon volume: 2ml. Sterile, for single use, in cartons of 10 units. | 5FG | R57422 | 2067145 |
| | 7FG | R57424 | 2067153 |
| HSG and HyCoSy Catheter with Stylet. HSG Catheter length 280mm, Maximum balloon volume: 2ml. Sterile, for single use, in cartons of 10 units. | 5FG | R57422-00-ST | 2067161 |
| | 7FG | R57424-00-ST | 2067170 |



CLIFFORDHALLAMHEALTHCARE

1300 720 274 www.ch2.net.au

Rocket Medical plc

Suite 6, 157 Gordon Street, Port Macquarie NSW 2444. AUSTRALIA

Tel: +61 (0) 431 146448

E:aus@rocketmedical.com

W:aus.rocketmedical.com

 Sales_007_ST AUS
 Rev: 11
 Date: 2021-09-09
 Copyright© 2014-21