




# Rocketmedical

## HSG & HyCoSy Catheters



Reproductive Medicine



**Rocket Medical HSG and HyCoSy Catheters have these well-established and trusted features:**

- **LATEX-FREE** balloon anchors the catheter and provides a contrast seal.
- **SMOOTH TIP AND LOW PROFILE BALLOON** ensures low trauma insertion and removal – ideal for the nulliparous cervix.
- **REMOVABLE STIFFENING SHEATH** to aid insertion in the heavily retroverted or anteverted cervix.
- **CONTRAST PATH TAP** ensures retention of contrast when further volumes of media are required.
- **INTEGRAL WIRE STYLET** for improved malleability without using tiresome removable stylets (ST versions)
  - ◇ Improves ease of insertion
- **COMPATIBLE** with aqueous based contrast media.

| Description  | Size | Order Code          | CH2 Code       |
|--|------|---------------------|----------------|
| <b>HSG and HyCoSy Catheter</b><br>HSG Catheter Length 280mm, Maximum balloon volume: 2ml.<br>Sterile, for single use, in cartons of 10 units.              | 5FG  | <b>R57422</b>       | <b>2067145</b> |
|  | 7FG  | <b>R57424</b>       | <b>2067153</b> |
| <b>HSG and HyCoSy Catheter with Stylet.</b><br>HSG Catheter length 280mm, Maximum balloon volume: 2ml.<br>Sterile, for single use, in cartons of 10 units. | 5FG  | <b>R57422-00-ST</b> | <b>2067161</b> |
|  | 7FG  | <b>R57424-00-ST</b> | <b>2067170</b> |

# CH2

CLIFFORDHALLAMHEALTHCARE

1300 720 274 [www.ch2.net.au](http://www.ch2.net.au)

**Rocket Medical plc**

Suite 6, 157 Gordon Street, Port Macquarie NSW 2444. AUSTRALIA

Tel: +61 (0) 431 146448

E:[aus@rocketmedical.com](mailto:aus@rocketmedical.com)

W:[aus.rocketmedical.com](http://aus.rocketmedical.com)

Sales\_007\_ST AUS  
Rev: 12  
Date: 2021-12-21  
Copyright© 2014-21