



CRAFT™ DUO-VAC™ illustrated

CRAFT™ Suction and CRAFT™ DUO-VAC™ Suction Pumps have been the gold standard for reliable oocyte aspiration for over 30 years. They have been continuously developed to provide smooth, low volume vacuum at pre-determined negative pressures, so allowing simple, low turbulence, oocyte recovery.

- **HIGHLY RELIABLE:** Rocket Medical CRAFT™ pumps have a proven 30 year history of reliability and excellent value.
- **SIMPLY, SAFE OPERATION** with air-operated foot switches on a solid moulded base. Lightweight connection set across the range to improve ease of handling and reliability.
- **LARGE, CLEAR GAUGE** with increased resolution improves setting accuracy.
- **ULTRA QUIET** diaphragm pump makes all Rocket Medical pumps easy to live with in the operating room.
- **FULLY VARIABLE VACUUM** from 0-200mmHg⁻¹. High Vacuum: 440mmHg⁻¹.



Description	Type	Order Code
CRAFT™ Suction Pump Complete with air operated vacuum footswitch, 5 single patient use water trap sets, mains supply cord, User Manual.	220-240VAC 50Hz	R29655
	110-120VAC 60Hz	R29654
CRAFT™ DUO-VAC™ Suction Pump Complete with air operated high/low vacuum footswitch, 5 single patient use water trap sets, mains supply cord, User Manual.	220-240VAC 50Hz	R29660
	110-120VAC 60Hz	R29661
Water Trap Set for CRAFT™ Suction Pumps 2.5m patient connection tube with male luer connector. Supplied sterile, for single use in cartons of 10 units.		R57685

