



CRAFT™ DUO-VAC™ illustrated

**CRAFT™ Suction and CRAFT™ DUO-VAC™ Suction Pumps** have been the gold standard for reliable oocyte aspiration for over 30 years. They have been continuously developed to provide smooth, low-volume vacuum at pre-determined negative pressures, so allowing simple, low-turbulence, oocyte recovery.

- **HIGHLY RELIABLE:** Rocket Medical CRAFT™ pumps have a proven 30 year history of reliability and excellent value.
- **SIMPLY, SAFE OPERATION** with air-operated foot switches on a solid moulded base. Lightweight connection set across the range to improve ease of handling and reliability.
- **LARGE, CLEAR GAUGE** with increased resolution to improve setting accuracy.
- **ULTRA QUIET** diaphragm pump makes all Rocket Medical pumps easy to live with in the operating room.
- **FULLY VARIABLE VACUUM**  
Medium Vacuum: 0-200mmHg<sup>-1</sup>. High Vacuum: 440mmHg<sup>-1</sup>.



Description	Type	Order Code
<b>CRAFT™ Suction Pump</b> Complete with air-operated vacuum footswitch, 5 single patient use water trap sets, mains supply cord, User Manual.	220-240V AC 50Hz	<b>R29655</b>
	110-120V AC 60Hz	<b>R29654</b>
<b>CRAFT™ DUO-VAC™ Suction Pump</b> Complete with air-operated high/medium vacuum footswitch, 5 single patient use water trap sets, mains supply cord, User Manual.	220-240V AC 50Hz	<b>R29660</b>
	110-120V AC 60Hz	<b>R29661</b>
<b>Water Trap Set for CRAFT™ Suction Pumps</b> 2.0m patient connection tube with male luer connector. Supplied sterile, for single use in cartons of 10 units.		<b>R57685</b>

