




- **TWO PART CONSTRUCTION IUI** with smooth outer sheath with high memory characteristics ensures ease of insertion, particularly in those patients where access to the external os is difficult.
- **SOFT INNER IUI** with fully domed distal tip and dual side eyes for optimum insemination.
- **COMPETITIVELY PRICED** to provide excellent value in routine IUI.
- **STERILE FOR SINGLE USE** - Packed in cartons of 25 units.



Rocket Medical DUO™ IUI Catheters are the latest generation of two part insemination catheters designed to improve ease of insertion and minimise uterine trauma.

Used throughout the world, Rocket Medical DUO™ Catheters are the device of choice when access is restricted and insemination likely to be difficult.

Specifications:

- (a) Inner Transfer Catheter: Polyurethane OD: 1.4 mm, ID: 0.9 mm, fully domed, dual side eye.
- (b) Outer Sheath: HDPE, 5 x 1 cm marking from distal tip.
- (c) Inner Catheter volume: 18cm: 0.5ml, 23cm: 0.63ml

Description	Length	Order Code
DUO™ IUI Catheters. Sterile, for single use, packed in cartons 25 units	18cm	R57625
	23cm	R57625-00-23