


Rocket Medical Embryon[®] Catheters

- **FIRM INNER CATHETER** for difficult transfers, with ultra smooth tip profile.
- **HIGH MEMORY OUTER SHEATH** for difficult cervical access.
- **REDESIGNED FINGER GRIP** gives unequalled grip and control whilst ensuring accurate placement of the catheter tip into the uterine cavity.
- **TOXICITY TESTING OF EVERY BATCH.** Rigorous quality controls ensure that all sensitive materials¹ are subject to LAL and independent single cell mouse embryo to blastocyst assay (MEA acceptance level: >80%) for the highest level of product safety.
- **STERILE FOR SINGLE USE** - Packed in cartons of 10 units.

References:

1. Craft I, Bernard A, Djahanbakan, Mcleod F (1982) Embryo transfer catheter material. Lancet. March 20, 1982, p68.

Rocket Medical Embryon[®] Embryo Transfer Catheters have provided reliable and secure transfer for over 25 years.

Used throughout the world, Embryon[®] Catheters are the device of choice when access is restricted and where transfer may be difficult.

Description	Length	Order Code
Embryon[®] ET Catheter (Peel Pouch) Sterile, carton 10 units	18cm	R57545
Embryon[®] ET Catheter (Tray Pack) Sterile, carton 10 units	18cm	R57570
Embryon[®] Stylet. Malleable, suitable for all new style Embryon [®] and Genesis [™] ET Catheters Sterile, carton 10 units	18cm	R57591-ET-18
	23cm	R57591-ET-23