

# Embryon<sup>®</sup> Genesis<sup>™</sup> Embryo Transfer



Embryon<sup>®</sup> Genesis<sup>™</sup> Catheter

- **FINGER GRIP DESIGN:** Retains our signature design for determining the depth of insertion of the inner catheter. The design provides unequalled grip and control whilst ensuring accurate placement of the catheter tip into the uterine cavity.
- **ULTRA-SOFT INNER CATHETER,** combined with an ultra-smooth tip profile, reduces the risk of junctional zone contraction.<sup>1</sup>
- **SMOOTH, THIN-WALLED OUTER SHEATH DESIGN** minimises cervical trauma. A new material, developed from established polymers, retains high memory characteristics despite its small diameter - ideal for patients where cervical access is difficult.
- **TOXICITY TESTING.** Rigorous quality controls ensure that all sensitive materials<sup>2</sup> are subject to LAL and independent single cell mouse embryo to blastocyst assay (MEA acceptance level: >80%). Batch tested to ensure the highest level of product safety.
- **STERILE FOR SINGLE USE** - Supplied in cartons of 10 units.

**References:**

1. Lensy P, Killick S.R et al (1998) Embryo transfer - can we learn anything new from observation of junctional zone contractions. Human Reproduction vol. 13, no.6, pp 1540-1546.
2. Craft I, Bernard A, Djahanbakan, Mcleod F (1982) Embryo transfer catheter material. Lancet. March 20, 1982, p68.

The Embryon<sup>®</sup> Genesis<sup>™</sup> catheter range has been developed using a soft inner catheter and ultra-smooth tip to help reduce the risk of junctional zone contractions which can displace or expel recently transferred embryos.



Description	Length	Order Code
<b>Embryon<sup>®</sup> Genesis<sup>™</sup> ET Catheter</b> Sterile, for single use, peel pouch supplied in cartons of 10 units	18cm	<b>R57630-00-18</b>
	23cm	<b>R57630-00-23</b>
<b>Embryon<sup>®</sup> Genesis<sup>™</sup> Trial Transfer Catheter</b> Sterile, for single use, peel pouch supplied in cartons of 10 units	18cm	<b>R57631-00-18</b>
	23cm	<b>R57631-00-23</b>
<b>Embryon<sup>®</sup> Stylet</b> Malleable, suitable for all new style Embryon <sup>®</sup> and Embryon <sup>®</sup> Genesis <sup>™</sup> ET Catheters Sterile, for single use, supplied in cartons of 10 units	18cm	<b>R57591-ET-18</b>
	23cm	<b>R57591-ET-23</b>