



- **COMPREHENSIVE KIT** containing the most commonly required items for fast, reliable bone marrow aspiration and biopsy.
- **REDUCE RISK AND STRESS** caused by searching for items in the ward or department, or by struggling with missing or worn instruments in reusable trays.
- **REDUCE WASTE:** combining components into a single kit significantly reduces wasted packaging.

**The Rocket Medical Bone Marrow Aspiration & Biopsy Kit** is the direct result of user research. The kit includes all items most commonly identified as required in the ideal ward or department procedure kit.

### KIT CONTENTS:

- 1 x T-handle bone marrow aspiration needle 15G
- 1 x trephine bone marrow biopsy Needle 11G
- 3 x 30mL Universal Containers, capped with labels
- 10 x VFM microscope slides, twinfrost, square edge
- 2 x 5-slot microscope slide transport boxes
- 3 x 2.5mL syringe, luer lock
- 1 x 10mL syringe, luer lock
- 1 x 20mL syringe, luer lock
- 1 x BD Microlance™ 3, 21g 40mm green needle
- 1 x Terumo Surguard 2, 25g 38mm orange safety needle
- 5 x 10cm x 10cm gauze swabs
- 1 x 60mL pot
- 1 x skin prep sponge stick
- 1 x 2-part fenestrated drape, 75 x 90cm
- 1 x square drape, 91 x 91cm<sup>2</sup>
- 1 x 5" sharp/blunt scissors
- 1 x Mepore self-adhesive wound dressing, 110mm x 150mm
- 1 x needle safe device (Atrion Needlewise)

### Ordering Information

### Code

#### Bone Marrow Aspiration and Biopsy Kit

Packed in procedure tray with sterile drape, in cartons of 5 kits.  
Sterile, for single patient use.

**R56750**