



Rocket® Single Lumen Aspiration Needles



I SX™ SPIRAL ECHO BEVEL DESIGN with high resolution echogenic surface preparation to give smoother, easier tissue penetration with extreme ultrasound definition.

I SHAPED FINGER GRIP with bevel marker and enhanced surfacing to significantly increase comfort and control.

I TWO NEEDLE STYLES:

I Rocket® Embryon®: dual angle bevel with 15mm brushed echogenic band

I Rocket® SX™: back sharpened bevel, 40mm high definition echo spiral.

I OXFORD TUBING SET: Single limb tube, set into a clear silicone bung with a 45° metal cannula connection (right).



Rocket® Spiral Echo Needles are the latest generation of high performance, high ultrasound definition oocyte aspiration sets from Rocket Medical plc.



| DESCRIPTION | DIMENSIONS | ORDER CODE |
|--|-------------------|--------------|
| Rocket® Embryon® Single Lumen Needles Sterile, in cartons of 10 sets | 17G x 32cm x 70cm | R57602 |
| | 16G x 33cm x 90cm | R57601-SX-90 |
| Rocket® SX™ Sterile, in cartons of 10 sets | 17G x 33cm x 70cm | R57602-00-SX |
| | 17G x 33cm x 90cm | R57602-SX-90 |
| | 17G x 35cm x 70cm | R57603-SX-70 |
| | 17G x 35cm x 90cm | R57603-SX-90 |